



CEREC Zirconia meso

- Pre-shaded blocks in classical colors
- Precise fit and correct occlusion thanks to enlarged milling and subsequent dense sintering
- Quick production as a result of the coordinated sintering programs in the CEREC SpeedFire
- The glaze firing gives a natural shine
- Fully anatomical zirconia guarantees high strength without the risk of chipping
- Kind to the gingiva thanks to the highest degree of biocompatibility and extraoral adhesion

Technical requirements

Production with: CEREC MC X, CEREC MC XL, CEREC MC XL Premium Package, inLab MC XL

Software: CEREC SW 4.5, inLab SW 16.0

Indication: Individually designed screw-retained crowns on Dentsply Sirona TiBase

Block	L × W × H	Blocks/Pack	TiBase Size	Ref. No.
CEREC Zirconia meso A1 S	24 × 23 × 21.5 mm	1	S	6623016
CEREC Zirconia meso A2 S	24 × 23 × 21.5 mm	1	S	6582428
CEREC Zirconia meso A3 S	24 × 23 × 21.5 mm	1	S	6582436
CEREC Zirconia meso A3.5 S	24 × 23 × 21.5 mm	1	S	6582444
CEREC Zirconia meso A1 L	24 × 23 × 21.5 mm	1	L	6623024
CEREC Zirconia meso A2 L	24 × 23 × 21.5 mm	1	L	6582451
CEREC Zirconia meso A3 L	24 × 23 × 21.5 mm	1	L	6583469
CEREC Zirconia meso A3.5 L	24 × 23 × 21.5 mm	1	L	6582477