

Our range of suction cannulas



Grey	Pack of 20	0700-059-00
Grey	Pack of 5	0700-059-50
Yellow	Pack of 5	0700-059-51
Blue	Pack of 5	0700-059-54

Against reflux: Protect universal cannula

Secondary air inlets can prevent the reflux of contaminated secretion from the suction system when suction is applied to the cannula. The special feature: Protect cannulas are particularly quiet, even with the inlets.



Grey	Pack of 4	0700-058-50
------	-----------	-------------

Prophylaxis cannula: with up to 2 x more aerosol absorption

Only from Dürr Dental: a specially developed rotating and fixable protective shield enables increased aerosol absorption and especially effective, ergonomic suction as part of preventive treatments with powder jet devices.



Grey	Pack of 20	0700-054-00
Grey	Pack of 5	0700-055-50
Yellow	Pack of 5	0700-055-51
Pink	Pack of 5	0700-055-53
Blue	Pack of 5	0700-055-54
Turquoise	Pack of 5	0700-055-55



For adults: universal cannula


Dürr Dental universal cannulas are particularly quiet. In a study by the Fraunhofer Institute, the measured volume was significantly below the level of comparable competitor products. They therefore ensure certified, low-noise suction of spray mist, and all of this without leaving behind any pressure sores.



Grey*	Pack of 5	A700056050
Yellow*	Pack of 5	A700056051
Pink	Pack of 5	A700056053
Blue*	Pack of 5	A700056054
Turquoise	Pack of 5	A700056055
ø 16 mm		

For children: Petito universal cannula

Children require a highly functional cannula that permits ergonomic treatment. The universal cannula Petito is smaller than the universal cannula; its Protect secondary air inlets provide optimal protection against reflux from the suction system. It has no corners or edges and leaves no sores.



ø 2.5 mm	Pack of 20	0700-007-50
ø 2.5 mm	Pack of 100	0700-007-51

Surgical suction cannula for single use

Surgical disposable cannulas with a 2.5 mm diameter and a cannula length of approx. 17 cm. The cannula has a tip with integrated Protect secondary air inlets. This offers the best possible protection from reflux from the suction system. The cannulas are individually packaged and sterile.