

## inCoris ZI meso

### For customized zirconium oxide abutments

The proven CAD/CAM production of individual zirconium oxide abutments offers you the opportunity to meet patient demands for natural tooth-colored implant-supported ceramic restorations.

- Zirconia blocks with prefabricated screw channel
- Two sizes and two shades
- After sintering, the mesostructure is adhesively connected to the titanium base (TiBase)
- Time savings compared to central production



inCoris ZI meso	
<b>Fabrication with</b>	CEREC MC X, CEREC MC XL, CEREC MC XL Premium Package, inLab MC XL
<b>Application</b>	Implant mesostructures and reduced crowns

inCoris ZI meso	REF.
inCoris ZI meso S F0.5	62 31 802
inCoris ZI meso L F0.5	62 31 810
inCoris ZI meso S F2	62 31 828
inCoris ZI meso L F2	62 31 836